REGENYAL  
**ARMONIOSO**  
CROSS-LINKED HYALURONIC ACID



1X3 ml

 **REGENYAL**  
LABORATORIES

REGENYAL  
**ARMONIOSO**  
CROSS-LINKED HYALURONIC ACID 1X3 ml

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The product is solely intended to be used by medical personnel. Do not use for applications other than indicated in this package leaflet.

#### DESCRIPTION

Regenyal Armonioso is a resorbable medical device (sterile, apyrogenic and physiological gel) to be used as filler for the correction of hypotrophy and to increase the volumes of the male and female genital areas. The main component is cross-linked hyaluronic acid of non-animal origin, produced by bacterial fermentation.

#### COMPOSITION

Cross-linked hyaluronic acid..... 75 mg/3 g  
Phosphate buffered saline, water for injectable solutions.....q.s.3 g

#### PACKAGE 1x3ml

- package leaflet
- sealed blister containing 1 single-dose/single-use pre-filled sterile syringe
- adhesive tabs to be applied to the patient card for product traceability.
- 1 needle and 1 cannula

#### METHOD OF ACTION

Regenyal Armonioso must be injected into the subcutaneous tissue for supplementation of the extracellular matrix and tissue volume increase.

#### INTENDED USE

Regenyal Armonioso is a medical device made in compliance with Directive 93/42/EEC as amended and added, intended for the treatment of male or female genital areas, in cases of hypotrophy and atrophy of the labia majora, volumetric increase of the pubis, volumetric increase of the penis. The outcomes of the procedure depend on the nature of the issues, results being the better the less noticeable the nature of the defect.

#### DIRECTIONS FOR USE

Prior to any procedure with Regenyal Armonioso, the doctor must conduct adequate anamnesis and an overall assessment of the patient's condition, to ensure the absolute absence of contraindications to the implant. The areas to be treated must be identified and assessed based on criteria to do with lines and symmetry one must adhere to.

The procedure may be performed under local anaesthesia in order to ensure the necessary comfort to the patient.

The physician should inform the patient beforehand on the methods of the procedure, its nature, warnings, precautions and possible individual outcomes, the potential adverse responses, the expected duration of treatment and the possible need for supplementary procedures for maintenance and/or detailed definition of the result achieved.

The area of the procedure should be cleaned with antiseptic solutions before injection.

Extract the syringe from the blister, remove the top as shown in the figure and screw the needle or the cannula to the luer-lock, complete with protection.

Only remove the cover just before performing the procedure.

Regenyal Armonioso is administered with a sterile needle-cannula compliant with Luer-Lock standards with standardised fittings, included in the package.

The injection should be performed in the superficial subcutaneous; the procedure is, however, at the physician's discretion and depends on the correction to be performed and the method implemented.

At the end of the procedure it is advisable to perform a gentle massage of the treated area for the product to be optimally distributed.

#### DOSE AND ADMINISTRATION

The volumetric markings on the syringes have indicative purposes, and the dosage to be used for the individual case is at the doctor's discretion. The frequency with which to repeat the procedure depends on several factors, concerning both the patient's physiology (skin type, individual metabolism, anatomy, age) and lifestyle. Another factor to be taken into account is related to the chosen implant techniques. The procedure should be repeated periodically in order to maintain the results achieved, every 5 - 6 months.

#### PRECAUTIONS

The usual precautions in case of percutaneous procedures must be taken during the procedure.

The risks are those of common infection related to the type of treatment. Regenyal Armonioso must not be used on patients suffering from:

- infectious or inflammatory processes close to the area of the procedure
- skin allergy in progress
- known hypersensitivity to Keloids
- allergy to components
- immune system disorders
- chronic pathological skin condition
- conditions affecting the clotting factors or in the case of anticoagulant therapy in progress.

Around the time of treatment, the patient should avoid taking substances (Aspirin, NSAIDs, Vit. E) that act on blood fluidity, in order to reduce to the minimum the likelihood of bruising or bleeding of the injected areas. Do not inject close to or during the menstrual cycle.

Use of Regenyal Armonioso should be carefully assessed by the physician on a case by case basis in patients who have been or are currently suffering from cancer in areas near the injection site and do not perform treatment in the case of progressive disease. The doctor must assure adequate surveillance of treated patients.

After the treatment avoid attending crowded and/or damp places (swimming pool, gym, sauna, Jacuzzi etc.) for at least 10 days.

Avoid sexual intercourse for 5days.

Avoid activities that cause pressure on the treated areas for a few days (cycling, horse riding etc.)

After use dispose of needles and syringes according to the indicated procedures for hospital waste.

#### ADVERSE REACTIONS, WARNINGS AND CONTRAINDICATIONS

As is the case for percutaneous injections, there may be phenomena of:

- inflammatory reaction (redness, swelling, etc.), sometimes associated with itching and pain at the injection site that usually resolve in a few days
- haematomas
- slight bleeding in treated areas
- hardening or lumps at the injection site
- skin colouring or discolouring at the injection site

These phenomena generally resolve in few days. Should they persist over a week or should other side effects appear, the patient must promptly contact their doctor.

Regenyal Armonioso must be used intracutaneously and must not be injected into blood vessels.

No overdosage phenomena or interaction with medicinal products are known.

Do not use during pregnancy or breastfeeding.

Regenyal Armonioso is available in sterile single use packages.

The product must be used immediately after opening.

It is forbidden to reuse the content for subsequent applications or on other patients.

Any unused product must be disposed of.

Do not use the product if the package is damaged.

Do not use jointly with other injectable products, or use on implants performed previously with other injectable products. Perform the treatment in a suitable facility following the appropriate techniques.

Fill in the adhesive labels included in the package, then apply one on the patient card stored at the doctor's surgery and hand the other over to the patient.

Keep away from the reach of children.

In rare cases, the product may cause an allergic reaction.

#### STORAGE METHODS

Regenyal Armonioso must be stored between 2°C and 28°C.

Do not freeze. Do not expose to heat sources.

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PIKDARE S.r.l
STERIMEDIX LTD
**EC REP** Bausch & Lomb GmbH,
**REGENYAL LABORATORIES Srl**